

Product datasheet for **GA215895**

Mouse Pkn3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pkn3 CRISPRa kit - CRISPR gene activation of mouse protein kinase N3
Symbol:	Pkn3
UniProt ID:	Q8K045
Synonyms:	AW209115; BC034126
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)