

Product datasheet for **GA215904**

Mouse Mill1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Mill1 CRISPRa kit - CRISPR gene activation of mouse MHC I like leukocyte 1
Symbol:	Mill1
Components:	GA215904G1, Mill1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215904G2, Mill1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215904G3, Mill1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8HWA5
Synonyms:	5530400I18Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)