

Product datasheet for **GA215910**

Mouse Sowahc activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Sowahc CRISPRa kit - CRISPR gene activation of mouse sosondowah ankyrin repeat domain family member C
Symbol:	Sowahc
Components:	GA215910G1, Sowahc gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215910G2, Sowahc gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215910G3, Sowahc gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8C0J6</u>
Synonyms:	4921515A04Rik; 9430099J10Rik; Ankrd57; C820004L04Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)