

Product datasheet for **GA215918**

Mouse Sh3pxd2b activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Sh3pxd2b CRISPRa kit - CRISPR gene activation of mouse SH3 and PX domains 2B
Symbol:	Sh3pxd2b
Components:	GA215918G1, Sh3pxd2b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215918G2, Sh3pxd2b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215918G3, Sh3pxd2b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A2AAY5
Synonyms:	fad49; G431001E03Rik; TSK4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)