

Product datasheet for **GA215931**

Mouse E030025P04Rik activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	E030025P04Rik CRISPRa kit - CRISPR gene activation of mouse RIKEN cDNA E030025P04 gene
Symbol:	E030025P04Rik
Components:	GA215931G1, E030025P04Rik gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215931G2, E030025P04Rik gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215931G3, E030025P04Rik gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	Gm664
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)