

Product datasheet for **GA215933**

Mouse Bahcc1 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Bahcc1 CRISPRa kit - CRISPR gene activation of mouse BAH domain and coiled-coil containing 1 |
| Symbol: | Bahcc1 |
| Components: | GA215933G1, Bahcc1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215933G2, Bahcc1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215933G3, Bahcc1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID: | Q3UHR0 |
| Synonyms: | B930044J06; BC060615 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online »](#)