

Product datasheet for **GA215935**

Mouse Sntg2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Sntg2 CRISPRa kit - CRISPR gene activation of mouse syntrophin, gamma 2
Symbol:	Sntg2
Components:	GA215935G1, Sntg2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215935G2, Sntg2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215935G3, Sntg2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q925E0</u>
Synonyms:	2210008K22Rik; 9530013L23; BB121248; G2syn; Syn5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)