

Product datasheet for **GA215938**

Mouse Tmem229b activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem229b CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 229B
Symbol:	Tmem229b
Components:	GA215938G1, Tmem229b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215938G2, Tmem229b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215938G3, Tmem229b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8BFQ2
Synonyms:	6330442E10Rik; 9830138K24; AU019070
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)