

Product datasheet for **GA215946**

Mouse Arhgef40 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Arhgef40 CRISPRa kit - CRISPR gene activation of mouse Rho guanine nucleotide exchange factor (GEF) 40
Symbol:	Arhgef40
Components:	GA215946G1, Arhgef40 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215946G2, Arhgef40 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215946G3, Arhgef40 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3UPH7</u>
Synonyms:	E130112L23Rik; Gm669; Solo
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)