

Product datasheet for **GA215948**

Mouse Carmil3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Carmil3 CRISPRa kit - CRISPR gene activation of mouse capping protein regulator and myosin 1 linker 3
Symbol:	Carmil3
Components:	GA215948G1, Carmil3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215948G2, Carmil3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215948G3, Carmil3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3UFQ8</u>
Synonyms:	Lrrc16b; mFLJ00240
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)