

## Product datasheet for **GA215961**

### Mouse *Atp13a5* activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                  |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                     |
| Description:  | Atp13a5 CRISPRa kit - CRISPR gene activation of mouse ATPase type 13A5                                                                                                                                                                                                                                           |
| Symbol:       | Atp13a5                                                                                                                                                                                                                                                                                                          |
| Components:   | GA215961G1, Atp13a5 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA215961G2, Atp13a5 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA215961G3, Atp13a5 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q3TYU2</a></u>                                                                                                                                                                                                                                                                                    |
| Synonyms:     | C630015F21Rik                                                                                                                                                                                                                                                                                                    |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                          |

#### Product images:



[View online »](#)