

Product datasheet for **GA215975**

Mouse Sap130 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Sap130 CRISPRa kit - CRISPR gene activation of mouse Sin3A associated protein
Symbol:	Sap130
Components:	GA215975G1, Sap130 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215975G2, Sap130 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215975G3, Sap130 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BIH0</u>
Synonyms:	2610304F09Rik; 6720406D06
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)