

Product datasheet for **GA216011**

Mouse Lrig2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lrig2 CRISPRa kit - CRISPR gene activation of mouse leucine-rich repeats and immunoglobulin-like domains 2
Symbol:	Lrig2
Components:	GA216011G1, Lrig2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216011G2, Lrig2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216011G3, Lrig2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q52KR2
Synonyms:	4632419I10Rik; BB096938; LIG-2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)