

Product datasheet for **GA216023**

Mouse Plekhg5 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Plekhg5 CRISPRa kit - CRISPR gene activation of mouse pleckstrin homology domain containing, family G (with RhoGef domain) member 5
Symbol:	Plekhg5
Components:	GA216023G1, Plekhg5 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216023G2, Plekhg5 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216023G3, Plekhg5 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q66T02
Synonyms:	BC023181; mKIAA0720; RP23-284C20.1; Syx
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)