

Product datasheet for **GA216025**

Mouse Pank4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pank4 CRISPRa kit - CRISPR gene activation of mouse pantothenate kinase 4
Symbol:	Pank4
UniProt ID:	Q80YV4
Synonyms:	D030031112Rik; R75150
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)