

Product datasheet for **GA216035**

Mouse Golga3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Golga3 CRISPRa kit - CRISPR gene activation of mouse golgi autoantigen, golgin subfamily a, 3
Symbol:	Golga3
UniProt ID:	P55937
Synonyms:	5330413L04; 5430416E01Rik; AI449376; AW490576; G1-499-14; Mea-2; Mea2; repro27
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)