


Product datasheet for **GA216047**

Mouse Clec4a1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Clec4a1 CRISPRa kit - CRISPR gene activation of mouse C-type lectin domain family 4, member a1
Symbol:	Clec4a1
Components:	GA216047G1, Clec4a1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216047G2, Clec4a1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216047G3, Clec4a1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	BC049354; Dcir4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)