

Product datasheet for **GA216060**

Mouse Ttl13 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ttl13 CRISPRa kit - CRISPR gene activation of mouse tubulin tyrosine ligase-like family, member 13
Symbol:	Ttl13
Components:	GA216060G1, Ttl13 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216060G2, Ttl13 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216060G3, Ttl13 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A4Q9F6</u>
Synonyms:	1700111A04Rik; 4921538O11; BC024578
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)