

Product datasheet for **GA216093**

Mouse **Taf1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Taf1 CRISPRa kit - CRISPR gene activation of mouse TATA-box binding protein associated factor 1
Symbol:	Taf1
UniProt ID:	<u>Q80UV9</u>
Synonyms:	250kDa; AU015687; B430306D02Rik; Ccg-1; Ccg1; KAT4; N-TAF1; Taf2a; TAFII250
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)