


Product datasheet for **GA216104**

Mouse Rp1l1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Rp1l1 CRISPRa kit - CRISPR gene activation of mouse retinitis pigmentosa 1 homolog like 1
Symbol:	Rp1l1
Components:	GA216104G1, Rp1l1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216104G2, Rp1l1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216104G3, Rp1l1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8CGM2</u>
Synonyms:	BC031365; BE983540; Dcdc4; Rp1hl1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)