


Product datasheet for GA216123

Mouse Poln activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Poln CRISPRa kit - CRISPR gene activation of mouse DNA polymerase N
Symbol:	Poln
Components:	GA216123G1, Poln gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216123G2, Poln gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216123G3, Poln gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q7TQ07</u>
Synonyms:	POL4P
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)