

Product datasheet for **GA216136**

Mouse Eif5a activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | Eif5a CRISPRa kit - CRISPR gene activation of mouse eukaryotic translation initiation factor 5A |
| Symbol: | Eif5a |
| Components: | GA216136G1, Eif5a gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216136G2, Eif5a gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216136G3, Eif5a gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID: | P63242 |
| Synonyms: | AA410058; D19Wsu54e; eIF-4D; eIF-5A; eIF-5A-1; eIF-5A1; Eif4d; Eif5a1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online »](#)