


Product datasheet for GA216147

Mouse Tm4sf19 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tm4sf19 CRISPRa kit - CRISPR gene activation of mouse transmembrane 4 L six family member 19
Symbol:	Tm4sf19
Components:	GA216147G1, Tm4sf19 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216147G2, Tm4sf19 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216147G3, Tm4sf19 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	EG277203
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)