

Product datasheet for **GA216155**

Mouse Trp53i11 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Trp53i11 CRISPRa kit - CRISPR gene activation of mouse transformation related protein 53 inducible protein 11
Symbol:	Trp53i11
Components:	GA216155G1, Trp53i11 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216155G2, Trp53i11 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216155G3, Trp53i11 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q4QQM4</u>
Synonyms:	Tp53i11
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)