

Product datasheet for **GA216163**

Mouse Cyp4a12a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cyp4a12a CRISPRa kit - CRISPR gene activation of mouse cytochrome P450, family 4, subfamily a, polypeptide 12a
Symbol:	Cyp4a12a
Components:	GA216163G1, Cyp4a12a gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216163G2, Cyp4a12a gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216163G3, Cyp4a12a gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q91WL5
Synonyms:	BC025936; Cyp4a12
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

