

## Product datasheet for **GA216164**

### Mouse Fam205c activation kit by CRISPRa

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                                  |
|----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                                     |
| <b>Description:</b>  | Fam205c CRISPRa kit – CRISPR gene activation of mouse family with sequence similarity 205, member C                                                                                                                                                                                                              |
| <b>Symbol:</b>       | Fam205c                                                                                                                                                                                                                                                                                                          |
| <b>Components:</b>   | GA216164G1, Fam205c gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074)<br>GA216164G2, Fam205c gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074)<br>GA216164G3, Fam205c gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074)<br>1 CRISPRa-Enhancer vector, SKU GEI00056<br>1 CRISPRa scramble vector, SKU GEI00077 |
| <b>UniProt ID:</b>   | <u><a href="#">Q80YD3</a></u>                                                                                                                                                                                                                                                                                    |
| <b>Synonyms:</b>     | 4931430D02; FAM205CP                                                                                                                                                                                                                                                                                             |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                          |

#### Product images:

