

Product datasheet for **GA216166**

Mouse C2cd3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	C2cd3 CRISPRa kit - CRISPR gene activation of mouse C2 calcium-dependent domain containing 3
Symbol:	C2cd3
Components:	GA216166G1, C2cd3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216166G2, C2cd3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216166G3, C2cd3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q52KB6
Synonyms:	AU020772
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)