

Product datasheet for GA216172

Mouse Tmtc2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmtc2 CRISPRa kit - CRISPR gene activation of mouse transmembrane and tetratricopeptide repeat containing 2
Symbol:	Tmtc2
Components:	GA216172G1, Tmtc2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074) GA216172G2, Tmtc2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074) GA216172G3, Tmtc2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074) 1 CRISPRa-Enhancer vector, SKU GEI00056 1 CRISPRa scramble vector, SKU GEI00077
UniProt ID:	<u>Q56A06</u>
Synonyms:	8430438D04Rik; D330034A10Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

