

Product datasheet for GA216174

Mouse Wfikkn2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Wfikkn2 CRISPRa kit - CRISPR gene activation of mouse WAP, follistatin/kazal, immunoglobulin, kunitz and netrin domain containing 2
Symbol:	Wfikkn2
Components:	GA216174G1, Wfikkn2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216174G2, Wfikkn2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216174G3, Wfikkn2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q7TQN3</u>
Synonyms:	2610304F08Rik; AY100450; Gaspl; WFIKKNRP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

