

## Product datasheet for **GA216214**

### Mouse H4c3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Hist1h4c CRISPRa kit - CRISPR gene activation of mouse histone cluster 1, H4c
Symbol:	H4c3
UniProt ID:	<a href="#">P62806</a>
Synonyms:	Hist1h4a; Hist1h4b; Hist1h4d; Hist1h4f; Hist1h4h; Hist1h4i; Hist1h4j; Hist1h4k; Hist1h4m
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)