

## Product datasheet for **GA216220**

### Mouse H4c18 activation kit by CRISPRa

#### Product data:

**Product Type:** CRISPRa Kits

**Description:** Hist1h4n CRISPRa kit - CRISPR gene activation of mouse histone cluster 1, H4n

**Symbol:** H4c18


**Components:** GA216220G1, H4c18 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)  
GA216220G2, H4c18 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)  
GA216220G3, H4c18 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)  
1 CRISPRa-Enhancer vector, SKU GE100056  
1 CRISPRa scramble vector, SKU GE100077

**UniProt ID:** [P62806](#)

**Synonyms:** Hist1h4a; Hist1h4b; Hist1h4c; Hist1h4d; Hist1h4f; Hist1h4h; Hist1h4i; Hist1h4j; Hist1h4k

**Format:** 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)