

Product datasheet for GA216232

Mouse H2ac10 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Hist1h2af CRISPRa kit - CRISPR gene activation of mouse histone cluster 1, H2af
Symbol:	H2ac10
Components:	GA216232G1, H2ac10 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216232G2, H2ac10 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216232G3, H2ac10 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8CGP5
Synonyms:	H2a-221
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

