

Product datasheet for **GA216233**

Mouse H2ac20 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Hist2h2ac CRISPRa kit - CRISPR gene activation of mouse histone cluster 2, H2ac
Symbol:	H2ac20
UniProt ID:	Q64523
Synonyms:	H2a-613b
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)