

Product datasheet for GA216235

Mouse H2bc3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Hist1h2bb CRISPRa kit - CRISPR gene activation of mouse histone cluster 1, H2bb
Symbol:	H2bc3
UniProt ID:	Q64475
Synonyms:	H2b-143
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)