

Product datasheet for **GA216273**

Mouse *Adgrl3* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Adgrl3</i> CRISPRa kit - CRISPR gene activation of mouse adhesion G protein-coupled receptor L3
Symbol:	<i>Adgrl3</i>
UniProt ID:	Q80TS3
Synonyms:	5430402I23Rik; C1RL-3; D130075K09Rik; Gm1379; LEC3; Lphn3; mKIAA0768
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)