

Product datasheet for **GA216284**

Mouse Ppm1h activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Ppm1h CRISPRa kit - CRISPR gene activation of mouse protein phosphatase 1H (PP2C domain containing) |
| Symbol: | Ppm1h |
| UniProt ID: | <u>Q3UYC0</u> |
| Synonyms: | A430075L18Rik; ARHCL1; C030002B11Rik |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)