

## Product datasheet for **GA216287**

### Mouse Lrtm1 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                            |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                               |
| Description:  | Lrtm1 CRISPRa kit - CRISPR gene activation of mouse leucine-rich repeats and transmembrane domains 1                                                                                                                                                                                                       |
| Symbol:       | Lrtm1                                                                                                                                                                                                                                                                                                      |
| Components:   | GA216287G1, Lrtm1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA216287G2, Lrtm1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA216287G3, Lrtm1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q8BXQ3</a>                                                                                                                                                                                                                                                                                     |
| Synonyms:     | A930016D02Rik                                                                                                                                                                                                                                                                                              |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                    |

#### Product images:



[View online »](#)