

Product datasheet for **GA216288**

Mouse Insyn1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Insyn1 CRISPRa kit - CRISPR gene activation of mouse inhibitory synaptic factor 1
Symbol:	Insyn1
Components:	GA216288G1, Insyn1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216288G2, Insyn1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216288G3, Insyn1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8CD60</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)