

Product datasheet for **GA216301**

Mouse *Idi1* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Idi1</i> CRISPRa kit - CRISPR gene activation of mouse isopentenyl-diphosphate delta isomerase
Symbol:	<i>Idi1</i>
UniProt ID:	P58044
Synonyms:	4832416K17Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)