

Product datasheet for **GA216316**

Mouse Fam168a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Fam168a CRISPRa kit - CRISPR gene activation of mouse family with sequence similarity 168, member A
Symbol:	Fam168a
Components:	GA216316G1, Fam168a gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216316G2, Fam168a gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216316G3, Fam168a gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BGZ2</u>
Synonyms:	2610030B18Rik; B930006L02Rik; mKIAA0280
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)