

Product datasheet for GA216327

Mouse Slc25a40 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Slc25a40 CRISPRa kit - CRISPR gene activation of mouse solute carrier family 25, member 40
Symbol:	Slc25a40
Components:	GA216327G1, Slc25a40 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216327G2, Slc25a40 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216327G3, Slc25a40 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8BGP6
Synonyms:	B230315F1IRik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

