

Product datasheet for **GA216330**

Mouse Cep295 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cep295 CRISPRa kit - CRISPR gene activation of mouse centrosomal protein 295
Symbol:	Cep295
Components:	GA216330G1, Cep295 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216330G2, Cep295 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216330G3, Cep295 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8BQ48
Synonyms:	5830418K08Rik; 5832426L23Rik; Gm1131; Kiaa1731
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)