

Product datasheet for GA216333

Mouse Ankar activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ankar CRISPRa kit – CRISPR gene activation of mouse ankyrin and armadillo repeat containing
Symbol:	Ankar
Components:	GA216333G1, Ankar gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074) GA216333G2, Ankar gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074) GA216333G3, Ankar gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074) 1 CRISPRa-Enhancer vector, SKU GEI00056 1 CRISPRa scramble vector, SKU GEI00077
UniProt ID:	<u>A2RT9I</u>
Synonyms:	4932422E22Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

