

Product datasheet for **GA216343**

Mouse Mpzl3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Mpzl3 CRISPRa kit - CRISPR gene activation of mouse myelin protein zero-like 3
Symbol:	Mpzl3
Components:	GA216343G1, Mpzl3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216343G2, Mpzl3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216343G3, Mpzl3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q3V3F6
Synonyms:	5430427F17Rik; A530065I17Rik; rc
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)