

Product datasheet for **GA216366**

Mouse Tmem229a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem229a CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 229A
Symbol:	Tmem229a
Components:	GA216366G1, Tmem229a gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216366G2, Tmem229a gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216366G3, Tmem229a gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	B9E1I9
Synonyms:	6332401O19Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)