

Product datasheet for **GA216368**

Mouse **Slc17a4** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Slc17a4 CRISPRa kit - CRISPR gene activation of mouse solute carrier family 17 (sodium phosphate), member 4
Symbol:	Slc17a4
Components:	GA216368G1, Slc17a4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216368G2, Slc17a4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216368G3, Slc17a4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5NCM1</u>
Synonyms:	9130214H05Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)