

Product datasheet for **GA216374**

Mouse Tmcc3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmcc3 CRISPRa kit - CRISPR gene activation of mouse transmembrane and coiled coil domains 3
Symbol:	Tmcc3
Components:	GA216374G1, Tmcc3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216374G2, Tmcc3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216374G3, Tmcc3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8R310</u>
Synonyms:	AW488095; C88213; C630016B22Rik; Tmcc1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)