

Product datasheet for **GA216377**

Mouse Dock6 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Dock6 CRISPRa kit - CRISPR gene activation of mouse dedicator of cytokinesis 6
Symbol:	Dock6
Components:	GA216377G1, Dock6 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216377G2, Dock6 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216377G3, Dock6 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8VDR9</u>
Synonyms:	2410095B20Rik; 4931431C02Rik; C330023D02Rik; mKIAA1395
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)