

Product datasheet for **GA216384**

Mouse Tns3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tns3 CRISPRa kit - CRISPR gene activation of mouse tensin 3
Symbol:	Tns3
Components:	GA216384G1, Tns3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216384G2, Tns3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216384G3, Tns3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q5SSZ5
Synonyms:	BC023928; F830010I22Rik; TEM6; Tens1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)