

Product datasheet for **GA216386**

Mouse A530016L24Rik activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	A530016L24Rik CRISPRa kit - CRISPR gene activation of mouse RIKEN cDNA A530016L24 gene
Symbol:	A530016L24Rik
Components:	GA216386G1, A530016L24Rik gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216386G2, A530016L24Rik gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216386G3, A530016L24Rik gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BNX7</u>
Synonyms:	Nrac
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)